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
Fig. 1 Fraction profile of Cr with InertSep® ME-1 250 mg/3 mL cartridge

Fraction test of InertSep® ME-1 was done independently for (a) Cr(III) and (b) Cr(VI) speciation. Break through: spiked 1 pm Cr(III) water (a) and spiked Cr(VI) water (b) loading into the separate SPE cartridges at pH 6; rinse: rinsed SPE cartridge with water at pH 6; 1st to 3rd elution solvent: 2 M HNO₃


Fig. 2 Fraction profile of Cr with MetaSEP AnaLig® Cr-02 200 mg/3 mL cartridge

Fraction test of MetaSEP AnaLig® was done independently for (a) Cr(III) and (b) Cr(VI) speciation. Break through: spiked 1 pm Cr(III) water (a) and spiked Cr(VI) water (b) loading into the separate SPE cartridges at pH 6; rinse: rinsed SPE cartridge with 0.1 M HCl; 1st to 3rd elution solvent: 4 M HNO₃

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
Fig. 3 Recovery of several metal ions in water by SPE using MetaSEP AnaLig® Cr-02

Each metal spiked at 200 ppb; Eluted by 4 M HNO₃; MetaSEP AnaLig® Cr-02 demonstrated by automated SPE workstation MetaPREP APS-1


Fig. 4 Evaluation of column flow rate for MetaSEP AnaLig® Cr-02 and AN-02 SPE cartridge prior to ICP-AES analysis

◆: MetaSEP AnaLig Cr-02 Cr(VI) recovery results; ●: MetaSEP AnaLig Cr-02 column breakthrough results; ▲: MetaSEP AnaLig AN-02 Cr(VI) recovery results; ■: MetaSEP AnaLig AN-02 column breakthrough results

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
Fig. 5 Sample preparation procedure for separation of Cr(VI) from Cr(III) matrix sample with highly selective resin, MetaSEP AnaLig® Cr-02, AN-02 SPE cartridge

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Table 1 Recovery test of Cr(VI) from actual sample with highly selective resin, SPE cartridges

Sample type	AnaLig® AN-02 Recovery, %	AnaLig® Cr-02 Recovery, %
A) Tap water	100.6	97.9
B) Plasma polishing waste water	91.5	96.9
C) Drill polishing waste water	89.2	78.1
D) Colored industrial wasted water	—	90.0
E) Colored soil extracts	—	88.5

(一) Not evaluated; Cr(VI) spiked level 100 ppb at A) sample; Cr(VI) spiked level 1 ppm for B), C), D), and E) samples; Detected by UV/VIS at 540 nm after MetaSEP AnaLig® AN-02, or Cr-02 SPE treatment

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
Fig. 6 Comparison of color from diphenyl carbazide reaction after MetaSEP AnaLig® Cr-02 treatment

① Industrial waste water sample, with non-SPE treatment; ② Industrial waste water sample after MetaSEP AnaLig Cr-02 treatment; ③ Cr(VI) standard solution after diphenyl carbazide reaction

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Determination of Cr(VI) by Solid-Phase Extraction with Highly Selective Resin

Yoshiaki FURUSHO¹, Masato ONO¹, Masayuki YAMADA², Takashi KITADE³ and Shoji MOTOMIZU⁴

¹ GL Sciences Inc., 6-22-1, Nishishinjuku, Shinjuku-ku, Tokyo 163-1130

² SII Nanotechnology Inc., RBM Tsukiji Bldg, 2-15-5, Shintomi, Chuo-ku, Tokyo 104-0041

³ M&S Instrument Inc., 2-12-4, Mikuni Honmachi, Yodogawa-ku, Osaka-shi, Osaka 532-0005

⁴ Chemistry and Biochemistry, Graduate School of Natural Science and Technology, Okayama University, 3-1-1, Tsushima-naka, Okayama-shi, Okayama 700-8530

(Received 23 December 2008, Accepted 19 January 2009)

Solid-phase extraction using molecular recognition technology (SPE-MRT) was done by using a chemically modified silica gel with covalently bound ligand polyatomic donor. Its highly selective adsorption behavior towards target elements was examined. The MetaSEP AnaLig[®] Cr-02 and AN-02 resins could adsorb Cr(VI), which contains chromate (CrO_4)²⁻ and dichromate (Cr_2O_7)²⁻ species, from a Cr(III) and Cr(VI) matrix solutions. Further investigation showed that the Cr(VI) species was quantitatively collected in the pH range of 0~6, whereas Cr(III) could not be retained on the resin in an acidic region (pH < 4). On the other hand, the InertSep[®] ME-1 synthesized methacrylate co-polymer resin with chelating function, imino-di-acetate group, could adsorb Cr(III) at a neutral region (pH 5~7), and did not adsorb selectively towards Cr(VI) at acidic to neutral region (pH 0~7). Determinations of Cr(III) and Cr(VI) were carried out by ICP-AES. Cr(VI) was also determined by using a diphenyl-carbazide colorimetric assay analysis with a UV/VIS spectrophotometer at 540 nm. To explore the merit of MRT resin, the speciation of chromium was demonstrated by using an automated SPE workstation, MetaPREP APS-1, whose metal parts were replaced with PEEK and a high-purity graphite tip needle. The SPE-MRT resins were packed in disposable polypropylene (PP) mini-columns and utilized as on-line connector devices with a sealing cap on the SPE cartridges, which were installed serially in the newly laboratory-assembled SPE automation. An Auto Pre-treatment System (APS) with ICP-AES and a UV/VIS spectrophotometer could provide highly sensitive determination for the speciation of Cr(III) and Cr(VI). The proposed system was successfully applied to the determination of Cr(VI) in spiked actual industrial wasted water and soil extract samples.

Keywords : solid-phase extraction ; SPE resin with molecular recognition technology ; chelating resin ; Cr(VI); ICP-AES ; diphenyl-carbazide colorimetric assay.